

Urine Drug Panels

	Panels:								
	17	16	15	14	12	10	9	7	5
Amphetamines amphetamine, MDA, MDEA, MDMA, methamphetamine	●	●	●	●	●	●	●	●	●
Cannabinoids carboxy-THC	●	●	●	●	●	●	●	●	●
Cocaine benzoylecgonine	●	●	●	●	●	●	●	●	●
Opiates codeine, hydrocodone, hydromorphone, morphine	●	●	●	●	●	●	●	●	●
Phencyclidine phencyclidine (PCP)	●	●	●	●	●	●	●	●	●
Benzodiazepines 7-aminoclonazepam, 7-aminoflunitrazepam, 2-hydroxyethylflurazepam, α -hydroxymidazolam, 7-aminonitrazepam, α -hydroxytriazolam, α -hydroxyalprazolam [‡] , nordiazepam [‡] , lorazepam [‡] , oxazepam [‡] , temazepam [‡]	●	●	●	●	●	●	● [‡]	● [‡]	
<small>‡Only denoted analytes included in 7- and 9-panels. All analytes included in 10-panel and above.</small>									
Barbiturates amobarbital, butalbital, pentobarbital, phenobarbital, secobarbital	●	●	●	●	●	●	●	●	
Methadone EDDP, methadone	●	●	●	●	●	●	●		
Propoxyphene norpropoxyphene	●	●	●	●	●	●	●		
Oxycodone oxycodone, oxymorphone	●	●	●	●	●	●			
Meperidine normeperidine	●	●	●	●	●				
Tramadol tramadol	●	●	●	●	●				
Fentanyl alfentanil, fentanyl, norfentanyl	●	●	●	●					
Sufentanil sufentanil (e.g. Sufenta [®]), norsufentanil	●	●	●	●					
Buprenorphine buprenorphine, norbuprenorphine	●	●	●						
Carisoprodol carisoprodol (e.g. Soma [®]), meprobamate (e.g. Equagesic [®])	●	●							
Zolpidem zolpidem (e.g. Ambien [®]), zolpidem phenyl-4-carboxylic acid	●								

Note: Ethanol is included in panels 10 and above.

Additional Add Ons and Available Options for Urine

Additional Add Ons and Available Options	Add On	Stand Alone
6-MAM 6-MAM	✓	X
Direct Ethanol Biomarkers ethanol, ethyl glucuronide, ethyl sulfate	✓	✓
Propofol Glucuronide propofol glucuronide (e.g. Diprivan®)	✓	✓
Kratom mitragynine	✓	✓



Urine collection materials. Collection instructions can be found at www.USDTL.com.

Urine provides the middle ground in drug testing, showing a history of drug exposure shorter than fingernails and hair, but longer than oral fluid (not offered at USDTL). A specimen of 10 ml provides information on the last 2-3 days of substance use history. We confirm all positive screening results by LC-MS/MS or GC-MS, offering urine testing panels with over 50 drug and alcohol biomarkers.